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mSoothe Foot massager

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MAX - F01DRS70A

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MAX - F01DRS70A

COMPANY PROFILE

After two decades of experience in the health industry, **maxcare** founder utilized his knowledge and expertise in massage therapy from an Asian-Pacific perspective to establish our brand - **maxcare**. Our massage therapy business was expanded to Hong Kong in 2012 and is conferred 2017 Hong Kong Top Brand Award. Based in Hong Kong, we are specializing in design & innovation, manufacturing, and retailing of massage and healthcare products.

Believing in the value of 'Technology improves our quality of life' and 'Caring for your family', **maxcare** has been innovating highly qualified products with international patents. Besides, **maxcare** products fulfill various international safety standards, namely EU CE Mark, US UL/ETL Mark, Germany GS Mark, TUV certification, and China CCC Mark.

Motivated by our value, **maxcare** founder has been leading us to be one of the most competitive brands in Hong Kong health industry. **maxcare** applies the contemporary fashion design – neat and modern – on the mechanical engineering of massage chairs. To deliver our value and offer incredible massage experience and healthcare solution to the public, **maxcare** devotes to the innovation, design, manufacturing the top-notch family massage products. By integrating the company philosophy – Spirit Of Healthy Life - into our products, **maxcare** brings health to everyone.

mSoothe Foot massager

Thank you for purchasing the Maxcare products. Before using the product, please read this manual thoroughly, particularly the Safety Precautions section, for the correct use of the product.

This manual should be well kept for easy reference. It enables you to perform a safe and efficient operation.

CONTENT

Safety Precautions	P.01 - P.03
Name and Function of Components	P.04
How to Use	P.05
Cleaning and Maintenance	P.06
FAQs and Solutions	P.07
Technical Specifications	P.08

Safety Precautions

1. Important safety warning

- Please use well-grounding power supply suited to this product.
- Please pull out the plug after use or before cleaning to avoid injury or damage to the product.
- Do not use unauthorized accessories.
- Do not use this product when the body is wet.
- Please read this manual instruction carefully before operation.
- Any other usage that is not listed in this manual instruction is forbidden.
- Please use the heating function with caution to avoid serious burnt.
- Recommended usage time of this product is 15 minutes/ session.
- Please do not keep this product running continuously for over 30 minutes.
- Please check if any part of this products is damaged before use. If so, please do not use it and contact us for repair.
- Please keep this product away from aerosol products or where oxygen is administered.
- Please do not remove the cover.
- Please do not use it if the cover is broken or damaged.
- Please do not use it if the shell or surface is damaged.
- Please do not allow foreign objects into this product.
- Please do not allow water or liquid into this product.
- Please do not use this product under a blanket or pillow. Overheating may occur and cause fire, electric shock or injury.
- Please do not cover this product when using it. The motor might overheat and cause malfunction or damages.
- Please do not fall asleep while using this product.
- Please do not use it if drunk or feeling unwell.
- Please do not use it after taking pain killers or sedatives.
- Please do not eat one hour before/ after using this product to avoid opposite effects.
- Please do not adjust to the maximum intensity at the beginning of the massage to avoid the risk of being injured.
- Keep the plug clean from needles, rubbish or water.
- Never cover the air opening when in use. Keep it clean from dirt, lint, hair and alike.
- If the power cord is damaged, do not attempt to change on your own. Please contact us for repair for safety reasons.

2. Environment for usage

- Please do not use it in a moist environment such as bathrooms.
- Please do not use it when there is a drastic change in temperature.
- Please do not use it under heavy-dusty or caustic environment.
- Please do not use it where there is not enough space or not good ventilation.
- Please switch off the power and unplug the machine after use.
- Please do not use it when surrounding temperature is up to 40°C.
- Please do not use it when surrounding temperature is below -20°C. The lubricant can thicken quickly and increase the friction between parts.
- Turn off and pull out the plug during thunderstorms.
- Do not use/leave this product outdoors.

Safety Precautions

- Do not place it near fire or directly under the sun.
- When moved from a low-temperature storage to a warm environment, the moisture inside the machine may affect, or damage, its functionality. Please rest the product within room temperature for 1 hour before use.

3. People who have the below symptoms, please consult a doctor before use

- People who have rashes or infectious skin diseases.
- People who have injuries or open wounds.
- People who have swollen or inflamed feet.
- People who are suffering from phlebitis, thrombosis, joint dysfunction, inflammation, weak bones or osteoporosis, etc.
- People who are suffering from heart disease or carrying electronic medical equipment such as heart pacemaker.
- People who have pins or screws inside the body, or artificial joints.
- People who have diabetes, sensory impairment, acute illnesses, heart diseases, or severe high blood pressure, etc.
- People who have a fever.
- People who have an abnormal or curved spine.
- Pregnant women or who has recently given birth.
- People who have injuries or body surface conditions.
- Children under 14-year old or people mentally unbalanced without being supervised are not allowed to use this machine.
- People who are told by doctor to have rest or people who feel unwell.
- This product has heating function. People who are insensitive to heat must be careful when using the appliance.
- People who are receiving medical treatment, especially those who are physically unfit.
- People who are at an increased risk of developing blood clots.
- People who aim for recovery.

4. Safety

- Please do not place heavy objects on this product.
- Please do not stand while using this product. Place and remove feet in and from this product while seated.
- Do not use other therapy devices while using this product.
- This product is not to be used for medical purposes. Should diagnosis or treatment is needed, please consult a health professional.
- Only one person at a time.
- Recommended usage time of this product is 15 minutes/ session.
- Check if the voltage is suitable to the specification of this product.
- Do not use wet hands to pull the plug.
- Do not allow liquid into this product, for the avoidance of electric shock or damage to this product.
- Do not pull or drag the power cord when unplugging.
- Do not damage the wires or change the circuit of this product.

Safety Precautions

- Do not alter the structure or parts of this product.
- Do not use wet cloths to clean the electric parts such as the switch and plug.
- Do not use any chemical solvent to clean this product. Use only a damp cloth to clean the surface.
- Stay away from this product during power cut to prevent injury in case power resumes suddenly.
- Stop using this product if there appears any irregularity, and please contact us immediately.
- Stop using this product if it causes discomfort, and please consult a doctor.
- This appliance is not intended for people (including children) with reduced physical, sensory or mental capabilities, or a lack of experience and knowledge, unless they are under supervision or instruction concerning the use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the appliance.
- If the power cord is damaged, it must be replaced by us in order to avoid a hazard.
- Do not carry this product by its power cord or use it as a handle.
- Keep the power cord away from hot or moving objects.
- Please cut off power and pull out the plug before moving.
- If the product is found damaged, broken, suffering from power leakage, or with exposed parts, please stop using it.

5. Maintenance

- This product should be maintained by us, users are prohibited to disassemble or maintain by themselves.
- Please always turn off the main power after use.
- Do not use this product if the socket is loose.
- If this product is to be left idle for long, please tie up the power cords and keep it in dry and dustless environment.
- Do not keep this product at a high temperature or near fire. Do not keep under direct sunlight.
- The power cord must be distanced from the heating surface.
- Please clean this product with dry cloths. Do not use thinner, benzene or alcohol.
- The mechanical components of this product are specially designed and made, no special maintenance required.
- Do not point sharp objects at this product.
- Please use it intermittently; do not keep it continuously running for long.

6. Solution of ordinary malfunctions

- It is normal if the motor makes noises during operation.
- If the rated working time is over, the power will automatically go off. If this product works continuously for too long, the temperature control will automatically turn off the machine, if so, please allow a 30-minute break before resuming.

Name and Function of Components

Product Structure

Touch panel ①

Touch panel ②

Removable and washable cover

Product features:

1. This product is developed by professional product R&D team. The design is elegant and with assorted colors. Classy lacquer surface radiates fashionable qualities. It is also light-weighted and portable.
2. Ergonomic design enables you to enjoy professional leg & foot massage at any time. The ergonomic design, unique styling and acupoint distribution all fit to the human leg & foot structure. Just like bringing a professional massagist home, you will enjoy ultimate relaxed massage on your leg & foot while watching TV program, or before sleep.
3. Double safety protection makes you feel secure when use this product. When the product is overheating or the current is overload, the circuit will break off automatically. This double protection is safe and reliable, ensure you to enjoy the massage joyous and carefree.
4. Combined various massage ways of kneading, scrapping, shiatsu, airbag press, warm foot etc.
5. Unique calf massage and care function.
6. Removable and washable cover, keep you clean at all time.
7. Equipped with double sensitive touch panel.
8. Use this massager 1-2 times a day, 5-15 minutes each time would effectively eliminate muscle fatigue, improve blood circulation. Use this massager for 5-15 minutes before sleep, would effectively improve sleep quality.

 HOW TO USE

Plug in the power, the machine into standby, ON / OFF button red light flashes.

Buttons	Description	Display
ON/OFF 	Press the ON/OFF button on the front panel of the massager to start the massager and enter the automatic working state.	Lights On
	Press for turn off	Light Twinkling
HEATING 	Press the HEATING button to turn on the foot heating function.	Lights On
	Press the HEATING button to turn off the foot heating function.	Light Off
STRENGTH 	Air pressure strength strong	Red Light
	Air pressure strength medium	Blue Light
	Air pressure strength weak	Purple Light
	Press the button for one second, you can turn off all the air pressure function.	Light Off
LEGS 	Foot rolling massage and air pressure work at the same time	Red Light
	Strengthen the leg massage, air pressure only	Blue Light
	Strengthen the foot scrapping function, leg air massage, the foot scrapping function still working	Purple Light
	When red light is on, press for one second to turn off massage function and only foot scrapping function	Light Off
SCRAPING 	Foot air pressure, roller work at the same time	Red Light
	Strengthen the foot scrapping function	Blue Light
	Strengthen the foot air pressure function	Purple Light
	Hold the button for one second, turn off the foot massage function	Light Off
VIBRATION 	Strong foot vibration function	Red Light
	Medium foot vibration function	Blue Light
	Soft foot vibration function	Purple Light
	Hold the button for one second, turn off the foot massage function	Light Off

Instructions:

- Place the massager on a clean floor and check if the machine's foot cloth cover is loose or broken
- Connect the adapter with the product, the other end connect to the power, the massager power indicator light starts flashing, the massager into the standby mode, people sitting in a chair and put the feet into the product(note not stand to use).
- Choose leg-and-foot massage mode or foot massage mode only by folding and unfolding the product. Folded product can massage the foot. Unfolded product can massage the foot and calf at the same time.
- Touch panel ① is invalid when the flip is folding, Touch panel ② is invalid when the flip is unfolding.
- If you want to change the flip state in the using, you need to press the ON/OFF button, until the product stops working and then lift the flip.
- In the working state press the power button or the end of the timer time, the massager will stop working, the power indicator flashes, enter the standby mode.

 **Cleaning and Maintenance****Cleaning**

- Before cleaning, please make sure the switch has been turned off, and the power plug has been pulled out from socket.
- Please use a piece of damp cloth to clean this product. The cloth can be soaked in soap water first, but please dry it to damp before cleaning.

**CAUTION**

- Do not pull the plug from socket with wet hands.
- Rinsing this product with water is not allowed.
- Do not use corrosive cleaning agents like benzene or thinner to clean the product, otherwise, the color of product will fade, and the internal parts may be corrupted or broken.
- Do not disassemble the product when it is wet.
- When not in use, turn off the main power, and remove the plug from socket.
- Keep the product clean and away from high temperature/ humid areas, or any environment with volatile gases.
- Keep the product in a carton, and put it in a safe, dry, well-ventilated and shaded place.
- If this product is to be left idle for long, properly cover it with a dust cover (self-purchase).
- Please keep the product away from direct sunlight, or high temperature, to avoid fading.

Kindly reminding

- If the product has been used for long, please check the conditions as below:
 - * Burnt smell.
 - * Unstable power supply after contact with the power cord.
 - * Power cord heating up.
 - * Power cords damage.
 - * Other abnormalities.

**CAUTION**

- If above mentioned happens, please stop using the product.
- Switch off the power supply, pull out the plug, and contact us.
- Do not disassemble the product by yourself.

FAQs and Solutions

No.	Issues	Possible Causes	Solution
1	"Tik tik" or "tak tak" sound during motion	Normal mechanical sound during motion or the frictional sound between the massage balls and cover	No action needed
2	Cannot be turned on	The adapter has not been well connected	Connect the adapter properly
		The power switch is off	Switch on the power button
		The power plug has not been well connected	Connect the plug to the socket properly
3	Abrupt stop during operation	Power disconnected	Check if the power cord is connected
		Rated working time (15 minutes) is up and the machine is automatically turned off	Press the power button to restart
		The machine has been in use for too long, causing overheating protection	Press the power button to restart the machine after cooling down
		Too much pressure on the machine and the massage balls are stuck	Relax the feet to lessen pressure on the machine

If the abnormal situations remain unsolved after taking above solutions, please contact us.

Technical Specifications

Product Name	: mSoothe Foot Massager
Product Model	: MAX-F01DRS70A
Rated Voltage	: 220V-50Hz
Rated Power	: 55W
Rated time	: 15 mins
Product Gross Weight	: 15 kg
Product Net Weight	: 11.5 kg
Product Size (Erected)	: 44 x 44 x 26 cm
Product Size (Stretched out)	: 44 x 44 x 45 cm
Execution Standard	: GB4706.1-2005 GB4706.10-2008

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Trademark holder: Maxcare Industrial Limited