

6In 1

加強版

脈衝按摩儀

6 In 1 Function Pulse Massager 2.0

MAX-M01MC0283
Instruction Manual



maxcare®



Note

- You can only adjust the mode and intensity level when the device is on, and connected to the 4 gel pads which are in contact with the skin. The device is programmed to turn off automatically after 30 seconds if the gel pads are not in contact with the skin.
- Each gel pad is good for about 100 uses. Do not place and remove from the skin if you do not intend to use it.
- If the gel pads are not placed securely on your skin, you will feel an uncomfortable tingling sensation. This may also happen if your skin is too oily or wet, or if the gel pads are worn and need to be replaced.


Care & Maintenance

- Do not store the device in humid or hot places
- Do not immerse the device or connection wire in water.
- If the gel pads are dirty or the adhesive is no longer sticky, wipe them gently using a damp cloth, or place them under gentle running water. Do not use harsh soaps. Allow to dry thoroughly before use.

Recommended Usage

- Recommended duration per area is 10-15 minutes.
- Do not use the massager for more than 30 minutes each time, 1 or 2 times per day.

Description & Function of Parts




How to Operate

1. Open the battery cover of the electronic device by sliding it gently downwards.
2. Inset two AAA battery into the battery compartment and slide the battery cover back into position.
3. Connect the connection wire to the output jack of the device.
4. Connect all 4 electrode gel pads to the snap button fasteners on the connection wire.
5. Remove the protective film of the gel pads and place the gel pads on the targeted body area.
6. Push the "ON/HI" button to turn on the device only when the gel pads are in position.
7. Press the "Mode" button repeatedly to select the desired programs. There are 6 available programs.
8. Press the "ON/HI" button slowly to increase the intensity, or the "OFF/LO" button to decrease the intensity. The selected intensity level will be displayed on the LCD.
9. Press "TIME" to set the timer for auto-off.
10. To switch the device off, press the "OFF/LO" button repeatedly until it turns off.
11. Replace the protective film back onto the gel pads after use. Remove the battery from the electronic device if it's not in use for an extended period of time.

Problem Shooting Guide

Problem	Possible Cause	Solutions
No display shown	<ul style="list-style-type: none">• Battery not installed• Battery installed incorrectly• Battery is dead	<ul style="list-style-type: none">• Install battery• Check the polarity of the battery• Replace battery
With display but not functioning	<ul style="list-style-type: none">• Disconnect with skin• Intensity level at 0• Microprocessor halted• Battery low	<ul style="list-style-type: none">• Place gel pad firmly onto skin• Push 'ON/HI' button to increase intensity• Reset by removing battery for 1 minute; reinsert battery• Replace battery
LCD appears normal but there is no pulsation	<ul style="list-style-type: none">• Intensity set to low level• Insufficient conductive gel on worn pad• The gel pads are positioned incorrectly	<ul style="list-style-type: none">• Push "ON/HI" button several times to increase intensity slowly• Clean gel pad or replace with new ones• Switch off the device. Move gel pads to other positions to locate trigger points



Product Specifications

Product Name:	6 In 1 Function Pulse Massager 2.0
Product Code:	MAX-M01MC0283
Voltage:	AAA battery x 2pcs
Frequency:	4-33 Hz
Programs:	6 modes, 10 intensity levels
Timer:	5/10/15/20/25/30 minutes
Output:	2.5 mA
Size:	10 x 7 x 25 cm
Weight:	55 g

Package includes:

- 6 in 1 Function Pulse Massager
- 4 Reusable electrode gel pads
- Connection wire with 4 snap button fasteners
- Instruction manual

Instruction Manual

- Please read this manual carefully before using 6 in 1 Function Pulse Massager.

- This is not a medical device. It is only intended for general physical conditioning. No medical claims are implied or warranted by the use of the product.

Caution

1. DO NOT use on the face, head, throat, chest or genital area.
2. DO NOT place on any wounds, or swollen, infected or inflamed areas.
3. DO NOT place on sensitive skin, or areas with skin conditions or nerve damage.
4. DO NOT use if pregnant, or on the abdomen during heavy menstrual flow, or immediately after childbirth.
5. DO NOT use after any recent surgical procedure, including Caesarian section.
6. DO NOT use while sleeping, driving or operating heavy machinery.
7. DO NOT use on children.
8. DO NOT use when there is a tendency to hemorrhage flowing acute trauma or fracture.
9. DO NOT use on areas with severe varicose veins.
10. DO NOT use if you suspect you have or are diagnosed with epilepsy, heart condition or multiple sclerosis, or if there is a risk of embolism.
11. DO NOT use if you have a pacemaker.